

Christopher Longmore¹, Nathalie Belot¹, Richard Clothier², Bushra Sim¹, Lottie Roscoe¹, Carol Treasure¹.
¹XCellR8 Ltd, The Innovation Centre, Sci-Tech Daresbury, Cheshire WA4 4FS, UK
²FRAME, Nottingham, UK
 Email: info@x-cellr8.com

INTRODUCTION

Acute toxicity is a key human health endpoint, requiring assessment by many global chemical safety regulators. Animal-based methods are used for acute toxicity with no validated alternative available or in development. XCellR8 has developed an acute toxicity screen using animal-product-free conditions. Comprised of a modified neutral red uptake (NRU) method and prediction model, we have screened 20 cosmetic ingredients to produce IC₅₀ values (concentration where cell viability is 50%). Through analysis of IC₅₀ values with the prediction model, 13/18 blind-tested method compatible chemicals were placed into the same GHS category as that derived from *in vivo* acute toxicity studies. This *in vitro* acute toxicity test is proposed as the initial step of an integrated approach to testing and assessment, incorporating multiple organ-specific endpoints. As a next step, XCellR8 aims to model metabolism in the current test and promote its utility for other chemical sectors.

METHODS

CELL CULTURE

Human dermal fibroblasts (HDFn) in Fibrolife™ Xeno-free medium – maintained at 37 °C/5% CO₂




Figure 1: Images of HDFn-XF cells at passage 1, day 1 (A), day 2 (B), and day 4 (C) x10 magnification

RESULTS

Table 1: IC₅₀ values of 20 test items and their *in vivo/vitro* toxicity categories. Items highlighted green showed concordance between classifications.

* Methylparaben and propylparaben had very close IC₅₀ and IC₂₀ values, so were classified as false positives. The closer these values are together, the more likely the potential for toxicity

Test Item	GHS oral tox. classification	<i>In vitro</i> acute tox. classification	IC ₅₀ Value (µg/ml)
Ammonium Lauryl Sulfate	4	4	129
Sodium Laureth Sulfate	4	4	150
Cocamidopropyl Betaine	4	4	67.8
Stearic Acid	5	4 or above (solubility limitation)	>100
Benzyl Acetate	5	5	>2000
Lemon Oil	5	4	306
Limonene	5	N/A (Too volatile)	<8000
Pulegone	4	5	>2000
Methylparaben	5	4*	352
Propylparaben	N/A	4*	155
Sodium Benzoate	5	5	22187
Phenoxyethanol	4	5	2205
Tudor Ebony (Black 11)	5	4 or higher	>100
Tudor Paeony (Red 7)	5	4 or higher	>200
Titanium Dioxide-Nano	N/A	5 or higher	>2000
Octyl Methoxycinnamate	5	5	>2000
Dimethicone	5	5	>200000
Glycerin	N/A	5	147172
Lanolin	N/A	N/A (solubility limitation)	>20
Triethanolamine	5	5	3364

ASSAY

- Solubility** tests followed OECD Test Guidance document 129 for 20 ingredients of varying known toxicity levels from five functional cosmetic groups
- Cells** seeded into 96-well plates and treated with 8 different concentrations of each test sample. A range-finding test was followed by a main experiment to determine IC₅₀ (concentration where cell viability is 50%).
- Neutral Red uptake** by the cells was determined by reading the absorbance at 540nm with a spectrophotometer. The data was used to determine the IC₅₀, allowing each chemical to be placed into a GHS category via a prediction model.

	Fatal in contact to the skin	Toxic in contact to the skin	Harmful in contact to the skin	Potentially harmful in contact to the skin
GHS	Categories 1 & 2	Category 3	Category 4	Category 5
IC ₅₀	ND	<10 µg/mL	10-1000 µg/mL	>1000 µg/mL

Figure 2: Prediction model for *in vitro* acute toxicity following determination of IC₅₀ values of 20 test items




Figure 3: Evaluation of cellular viability by neutral red uptake for four test items

DISCUSSION

The results show 13 chemicals for which *in vitro* IC₅₀ values were determined and were placed into the same GHS classification as that based on LD₅₀ data obtained by oral acute toxicity. Removing the chemicals for which no IC₅₀ could be determined due to solubility limitations, this resulted in the test correctly classifying 13/18 chemicals (72%). Cytotoxicity is a key event in modes of action associated with acute health effects, covering mechanisms of toxicity common to most cell types that can lead to organ failure, including disruption of membrane function – the endpoint measured in the neutral red uptake test. We envisage that a robust human cytotoxicity assay would constitute an early-stage screen for highly toxic chemicals. Chemicals that are not classified as highly toxic in the cytotoxicity pre-screen would then be tested further using organ-specific endpoints. XCellR8 will seek to incorporate the study of metabolism into the assay, in addition to refining the prediction model. It is hoped this acute toxicity screen will form the first stage of an integrated *in vitro* testing strategy for acute toxicity.

REFERENCES

- Hoffmann *et al.* (2010). Acute oral toxicity: variability, reliability, relevance and interspecies comparison of rodent LD₅₀ data from literature surveyed for the ACuteTox project. *Regul Toxicol Pharmacol.* 58(3): 395-407.
- Kinsner-Ovaskainen, A. *et al.* (2013). Selection of test methods to be included in a testing strategy to predict acute oral toxicity: an approach based on statistical analysis of data collected in phase 1 of the ACuteTox project. *Toxicology In Vitro.* 27(4): 1377-94.