

# Short Communication

## Adaptation of the KeratinoSens™ skin sensitisation test to animal-product-free cell culture

Nathalie Belot, Bushra Sim, Christopher Longmore, Lottie Roscoe and Carol Treasure  
XCellR8 Ltd, the Innovation Centre, Sci-Tech Daresbury, Keckwick Lane, Daresbury, Cheshire, UK

### Summary

Skin sensitisation is the process by which a substance leads to an allergic reaction following skin contact. The process has been described as an adverse outcome pathway (AOP), including several key events, from skin penetration and covalent protein binding, to keratinocyte activation, dendritic cell activation and T-lymphocyte proliferation. The *in vitro* assay KeratinoSens™ measures the activation of keratinocytes. It is fully accepted at a regulatory level (OECD TG 442d) and appropriate for compliance with a range of legislation including the EU Cosmetics Regulation, REACH, and the CLP Regulation.

Currently, many *in vitro* methods use animal-derived components in the cell culture systems. Many stakeholders in the cosmetics industry have both scientific and ethical concerns relating to this issue and have stated a strong preference for fully human *in vitro* test systems. We have adapted the KeratinoSens™ method to animal product-free conditions, and carried out an in-house validation with 21 reference substances, including those listed in the Performance Standards associated with OECD TG442d. The modified method was shown to be totally equivalent to the Validated Reference Method (VRM), with comparable values for accuracy (85.7%), sensitivity (84.6%) and specificity (87.5%), and all acceptance criteria being met. In Europe, data generated by the adapted method may be used in REACH submissions, and we are now seeking approval to list the adaptation in OECD TG 442d, enabling formal compliance with a range of global regulations.

Keywords: skin sensitization, KeratinoSens™, alternatives to animal testing, 3R, *in vitro*

### 1 Introduction

Skin sensitization is a multi-step process which has been described as an Adverse Outcome Pathway (AOP). Following skin penetration, chemicals may bind covalently to skin proteins (Key Event 1), followed by the activation of the keratinocytes (Key Event 2) and dendritic cell activation (Key Event 3). Dendritic cells migrate to the lymph nodes and induce the proliferation of T-lymphocytes (Key Event 4). Under pressure from legislation, including the European Cosmetics Regulation 1223/2009 and REACH (Registration, Evaluation, Authorization and restriction of Chemicals), many groups have worked in recent years on identifying relevant alternatives to traditional animal methods to determine the skin sensitization potential of chemicals (Reisinger et al., 2015). Alternatives for Key Events 1-3 have recently been adopted as OECD guidelines: DPRA (OECD TG442c); KeratinoSens™ (OECD TG442d) and h-CLAT (OECD TG442e).

The KeratinoSens™ method addresses the second key event of skin sensitization - the activation of keratinocytes. KeratinoSens™ cells are derived from the human keratinocyte cell line, HaCat, containing a luciferase gene that is under the control of a constitutive promoter fused with the Antioxidant Response Element (ARE) from a gene that is known to be up-regulated by contact sensitizers. The majority of skin sensitizers induce this pathway and, therefore, the luciferase signal reflects the activation by sensitizers of endogenous Nrf2 dependent genes. This method was fully validated and gained regulatory acceptance in 2015 (Natsch et al., 2010; Natsch et al., 2011; Natsch et al., 2013).

The KeratinoSens™ method represents a great advancement in the replacement of animal testing. However, the published protocol currently includes the use of animal-derived cell culture components (bovine serum and porcine trypsin). XCellR8's mission is the full replacement of animal testing for cosmetics and their ingredients globally, motivated by an industry requirement to maximise the human relevance of *in vitro* models, along with animal welfare considerations. This strategy includes eradication of the significant number of animal-derived components that are still widely used in cell culture methods (including serum, tissue extracts and antibodies). In this project, minor adaptations to an existing method enable the complete replacement of animals for testing key event 2 in the skin sensitisation Adverse Outcome Pathway. They enable the method described in OECD TG 442d to be conducted without the need for a number of animal components by replacing them with human equivalents. Key animal components replaced are foetal calf serum and porcine trypsin. The significant animal welfare concerns around the production of these components has been well documented (Jochemset al., 2002).

Received January 31, 2017;  
Accepted March 16, 2017;  
Epub March 16, 2017  
<https://doi.org/10.14573/altex.1701311>



This is an Open Access article distributed under the terms of the Creative Commons Attribution 4.0 International license (<http://creativecommons.org/licenses/by/4.0/>), which permits unrestricted use, distribution and reproduction in any medium, provided the original work is appropriately cited.

Here, we describe the internal validation of the modified method using 21 reference chemicals, including those described in the Performance Standards (OECD, 2015) for skin sensitization and in OECD guideline Test Guideline 442d. The reference chemicals are comprised of a panel of different sensitizer categories (minimal to severe and non-sensitizers), as defined by the animal-based Local Lymph Node Assay (LLNA).

## 2 Materials and methods

### *Test Chemicals*

All chemicals were purchased from Sigma-Aldrich (UK) and selected following description in the Performance Standards (OECD, 2015) for skin sensitization and in OECD Test Guideline TG442d. Cinnamic aldehyde was used as the positive control in the assays (Sigma-Aldrich, UK CAS # 14371-10-9). Test chemicals were dissolved in dimethyl sulfoxide (DMSO) at a concentration of 200mM. They were further diluted in culture medium to a final concentration of 2000µM in 1% DMSO and then serially diluted in culture medium with 1% DMSO (DB-ALM Protocol 155 (DB-ALM, 2013)). The reference chemicals were: isopropanol (67-63-0); salicylic acid (69-72-7); lactic acid (50-21-5); glycerol (56-81-5); 4-methoxyacetophenone (100-06-1) chlorobenzene (108-90-7); methyl salicylate (119-36-8); sulfanilamide (63-74-1); ethylene glycol dimethacrylate (97-90-5); phenyl benzoate (93-99-2); eugenol (97-53-0); 2-mercaptobenzothiazole (149-30-4); citral (5392-40-5); isoeugenol (97-54-1); methyl dibromoglutaronitrile (35691-65-7); 4-methylaminophenol sulphate (55-55-0); parphenylenediamine (106-50-3); 2,4-dinitrochlorobenzene (97-00-7); 4-nitrobenzylbromide (100-11-8); oxazolone (15646-46-5); cinnamyl alcohol (104-54-1).

### *Cell culture*

The KeratinoSens<sup>TM</sup> cell line contains a stable insertion of a luciferase gene under the transcriptional control of a constitutive promoter fused with the ARE (Antioxidant Response Element) from the AKR1C2 gene (Lou et al., 2006). KeratinoSens<sup>TM</sup> cells were developed by Emter et al. (2010) and purchased from Givaudan. Upon first thawing from cryopreservation, cells were adapted to cell culture medium containing human serum. Cells were routinely maintained in Dulbecco's Modified Eagles Medium (DMEM) containing Glutamax (Life Technologies), 10% human serum (pooled human male AB plasma, 60-70 donors, Sigma-Aldrich UK) and Geneticin (500µg/ml, Life Technologies) at 37°C, 5% CO<sub>2</sub> to a maximum passage number of 22. They were grown to a maximum of 80-90% confluence. Sub-culture was performed using TrypZean® (Sigma-Aldrich) - a recombinant, animal-free equivalent of porcine trypsin.

### *KeratinoSens<sup>TM</sup> assay*

The KeratinoSens<sup>TM</sup> assay was performed as previously described (Emter, et al., 2010; DB-ALM, 2013), with the following adaptations and optimization under animal-product-free conditions.

Cells were seeded in 3 parallel replicate 96-well plates for luciferase activity testing and in one 96 well plate for MTT viability testing in medium containing 10% human serum. Plates were incubated at 37°C, 5% CO<sub>2</sub> for 24 hours prior to treatment. At the dosing time point, culture medium was removed and replaced with fresh culture medium containing 1% human serum and 1% DMSO. Cells were then exposed to the reference chemicals and controls, at 12 different concentrations (max 2000µM) and incubated for 48 hours at 37°C, 5% CO<sub>2</sub>. After the incubation period, luciferase activity was evaluated by luminescence measurement (Luciferase Assay System, Promega, UK) and cell viability using the MTT viability assay (Sigma-Aldrich). The assay was repeated 3 times. For each chemical in each repetition and at each concentration, the luciferase gene induction was compared to the negative control (1% DMSO). A maximum of 3 outliers out of 9 values were removed to obtain a CV value below 20%. As per the VRM (Validated Reference Method), the following parameters were calculated: the I<sub>max</sub> value (maximal average fold induction of luciferase activity observed at any concentration of the tested chemical and positive control); the EC<sub>1.5</sub> value, representing the lowest concentration for which induction of luciferase activity is above the 1.5 fold threshold (i.e. 50% enhanced luciferase activity); and the IC<sub>50</sub> and IC<sub>30</sub> concentration values for 50% and 30% reduction of cellular viability respectively.

The positive control, cinnamic aldehyde, was used at a concentration range of 8 to 128µM and required to meet the following acceptance criteria: average induction (I<sub>max</sub>) in the three replicates for cinnamic aldehyde at 32 µM between 1.6 and 3; EC<sub>1.5</sub> value between 16 µM and 32 µM.

Test items were considered as positive for skin sensitisation if the following conditions were met in 2 of 3 repetitions: The I<sub>max</sub> was higher than 1.5 fold and statistically significantly different as compared to the solvent (negative) control (as determined by a two-tailed, unpaired Student's T-test); the cellular viability was higher than 70% at the lowest concentration with induction of luciferase activity above 1.5 fold (i.e. at the EC<sub>1.5</sub> determining concentration); the EC<sub>1.5</sub> value was less than 1000µM; there was an apparent overall dose-response for luciferase induction.

## 3 Results

The KeratinoSens<sup>TM</sup> protocol adapted to animal-free conditions was performed to predict the skin sensitization potential of 21 reference chemicals, including those listed in the Performance Standards (OECD, 2015) and OECD Test Guideline 442d.

For each chemical, the I<sub>max</sub> value and EC<sub>1.5</sub> concentration were calculated (Tab.1) (I<sub>max</sub>, maximal average fold induction of luciferase activity observed at any concentration of the tested chemical and positive control); the EC<sub>1.5</sub> value, representing the lowest concentration for which induction of luciferase activity is above the 1.5 fold threshold (i.e. 50% enhanced luciferase activity).

**Tab. 1: Results obtained at XCELLR8 for the 21 reference chemicals (mean of 3 independent repetitions in triplicate)**

These results are compared to VRM KeratinoSens™ assay results and to *in vivo* classification of the chemicals, based on EC3 values from literature results of the LLNA test method: >10 weak; <10 and >1 moderate; <1 and >0.1 strong; <0.1 extreme (OECD N°213, 2015). The 3 underlined chemicals were misclassified as compared to *in vivo* results. Data from Emteret al. 2010 (1) and Natsch and Emter, 2008 (2), lmax= maximal fold gene induction; EC1.5: concentration in µM for 1.5-fold induction, n.i.: no induction.

Reference Substances	<i>In vivo</i> category (according to LLNA)	VRM prediction			XCELLR8 prediction		
		lmax (fold)	EC1.5 (µM)	Prediction	lmax (fold)	EC1.5 (µM)	Prediction
Isopropanol	Non Sensitizer	1.2 (1)	n.i. (1)	Non Sensitizer (1)	1.2	n.i.	Non Sensitizer
Salicylicacid	Non Sensitizer	1.1 (1)	n.i. (1)	Non Sensitizer (1)	1.4	n.i.	Non Sensitizer
Lactic Acid	Non Sensitizer	1.3 (1)	n.i. (1)	Non Sensitizer (1)	1.3	n.i.	Non Sensitizer
Glycerol	Non Sensitizer	1.2 (1)	n.i. (1)	Non Sensitizer (1)	1.4	n.i.	Non Sensitizer
<u>4-methoxy-acetophenone</u>	Non-Sensitiser	1.7(2)	449.3(2)	Sensitizer(2)	2.1	620	Sensitizer
Chlorobenzene	Non-Sensitiser	1.2(1)	n.i.(1)	Non-Sensitizer(1)	1.2	n.i.	Non-Sensitizer
Methyl salicylate	Non-Sensitiser	1.2(1)	n.i.(1)	Non-Sensitizer(1)	1.2	n.i.	Non-Sensitizer
Sulfanilamide	Non-Sensitiser	1.4(1)	n.i.(1)	Non-Sensitizer(1)	1.1	n.i.	Non-Sensitizer
Cinnamyl Alcohol	Sensitizer (Weak)	1.7(1)	123.6(1)	Sensitizer(1)	4.2	20	Sensitizer
Ethylene glycol dimethacrylate	Sensitizer (Weak)	188.4(1)	57.4(1)	Sensitizer(1)	4.8	29	Sensitizer
<u>Phenyl benzoate</u>	Sensitiser (Weak)	1.3(1)	n.i. (1)	Non-sensitizer(1)	1.1	n.i.	Non-Sensitizer
<u>Eugenol</u>	Sensitiser (Weak)	1.3(1)	n.i.(1)	Non-Sensitizer(1)	2.2	286	Non-Sensitizer borderline
2-Mercaptobenzothiazole	Sensitizer (Moderate)	8.8 (1)	48.1(1)	Sensitizer(1)	6.9	57	Sensitizer
Citral	Sensitiser (Moderate)	96.4(1)	23.2(1)	Sensitizer(1)	3.8	18	Sensitizer
Isoeugenol	Sensitiser (Moderate)	6.4(1)	16.1(1)	Sensitizer(1)	3.4	20	Sensitizer
Methyldibromoglutaronitrile	Sensitizer (Strong)	4(1)	7.8(1)	Sensitizer(1)	2.7	8	Sensitizer
4-Methylaminophenol sulphate	Sensitizer (Strong)	5.9(1)	9.4(1)	Sensitizer(1)	36.1	4	Sensitizer
Para-phenylenediamine	Sensitiser (Strong/Extreme)	26.8(1)	5(1)	Sensitizer(1)	28.2	6	Sensitizer
2,4 Dinitro-chorobenzene	Sensitizer (Extreme)	14.8(1)	2.5(1)	Sensitizer(1)	8.5	1	Sensitizer
4-Nitrobenzylbromide	Sensitiser (Extreme)	6.9(1)	1.3(1)	Sensitizer(1)	10.5	<0.98	Sensitizer
Oxazolone	Sensitiser (Extreme)	2.4(1)	175.5(1)	Sensitizer(1)	5.4	129	Sensitizer

Among the 8 reference compounds classified as non-sensitizers *in vivo* (based on published LLNA results), one was misclassified (4-methoxy-acetophenone), in common with the KeratinoSens™VRM.

Regarding the 13 chemicals classified as sensitizers *in vivo* (based on published LLNA results), 2 were misclassified (phenyl benzoate and eugenol), in common with the KeratinoSens™VRM.

Results show that the prediction for the 21 reference chemicals was identical to the VRM (Emteret al., 2010; Natsch and Emter, 2008).

Both assays correctly predicted 18 chemicals of the total 21 chemicals and incorrectly predicted the same 3 chemicals, when compared with LLNA results:4-methoxy-acetophenone was predicted as a skin sensitizer (“false positive”) while phenyl benzoate and eugenol were predicted to be non-sensitizers (“false negatives”). In compliance with the acceptance criteria of the Performance Standards, no strong or extreme sensitizer was under-predicted with the adapted method.

Furthermore, in the case of eugenol, the result obtained in our laboratory was a “borderline” classification, as the 2 last concentrations were determined as positive but with a cell viability below the 70% threshold.

## 4 Discussion

Many advances have been made, in terms of achieving full regulatory acceptance for scientifically advanced *in vitro* methods to replace traditional animal-based safety testing. Examples include: OECD TG431, OECD TG439 and OECD TG492 using reconstructed human tissue models to assess skin corrosion, skin irritation and eye irritation respectively. Recently, 3 methods for the assessment of key events in the skin sensitization adverse outcome pathway (AOP) have gained full acceptance and have been incorporated into the REACH Regulation in Europe as the methods of choice for this endpoint. These are OECD TG442c (Direct Peptide Reactivity Assay “*in chemico*” method for Key Event 1 – measuring protein binding); OECD TG442d (ARE-Nrf-2 Luciferase methods, including KeratinoSens™, for Key Event 2 – measuring keratinocyte activation); and OECD TG442e (h-CLAT, for Key Event 3 – measuring dendritic cell activation). The traditional animal-based method for the assessment of skin sensitisation, OECD TG429 - the Local Lymph Node Assay (LLNA) - may now only be used as a last resort, representing a significant shift towards the adoption of *in vitro* technologies in the context of regulatory safety testing.

In spite of this substantial progress, it is important to recognize that many *in vitro* methods currently use animal-derived components such as bovine serum, mouse antibodies, rat tissue extracts, and a variety of animal-derived cell lines. As such, these methods still ultimately require the sacrifice of animals and cannot be considered completely animal-free. Such an approach presents scientific limitations, by unnecessarily compromising the direct relevance of the test system to humans. The significant animal welfare concerns around the production of these components has also been well documented (Jochemset al., 2002; Even, et al., 2006). Importantly, these scientific and ethical drawbacks give rise to concerns within the cosmetics and personal care industry. The European Cosmetics Regulation 1223/2009 includes an animal testing ban for cosmetic products and ingredients. Globally, animal testing is a highly emotive issue in many countries, and there is increasing demand from consumers for “cruelty-free” cosmetic products, in which neither the formulation nor individual ingredients have been tested on animals. Many ethical cosmetic companies seek to go beyond the “bare minimum” expectation, and to satisfy consumer demand, by avoiding the use of animal-derived components altogether. In response to this demand, new methods developed at XCellR8 are animal-product-free from the outset, and the organization openly encourages other laboratories to adopt the same approach, avoiding the need for adaptation *after* the formal process of validation and regulatory acceptance has been completed. In addition, we have been asked by leading cosmetic companies to adapt existing regulatory methods to animal-product-free conditions, and to seek regulatory acceptance for the updated protocols.

In this project, the widely used KeratinoSens™ test method for skin sensitisation (OECD TG 442d) was adapted to animal-product-free conditions and further optimised by making the following minor adaptations to the established protocol: 1. The use of human serum in cell culture (OECD TG 442d states only “serum” but the Validated Reference Method (VRM) protocol states “foetal calf serum”). 2. Non-animal recombinant trypsin (TrypZean®) used to harvest cells. 3. MTT method adapted to the widely-used protocol consisting of 3 hours’ incubation, solubilisation in isopropanol and absorbance measurement at 570nm. 4. Positive control (cinnamic aldehyde) concentration range optimised to 8-128µM for the adapted test conditions.

The KeratinoSens™ cell line adapted well to routine culture in human serum (in place of bovine serum), showing healthy morphology and comparable growth rates. This enabled the rapid creation of an internal cell bank for long term use in the test. A maximum passage number of 22 was set, to reflect the optimum passage window for assay performance using the cells cultured in human serum.

Since there have been no changes to the cell lines or key steps in the protocols, the changes constitute an adaptation to OECD TG442d, and *not* a “me-too” method.

Subsequently, the adapted method was validated in-house using the 10 Proficiency Chemicals listed in OECD TG442d, and the 11 additional chemicals listed in the associated Performance Standards. (Of the 20 chemicals listed in the Performance Standards, 9 were already present in the Proficiency Chemicals listed in OECD guideline TG442d). These chemicals spanned the full range of levels of skin sensitisation observed in the commonly used animal-based test for skin sensitisation (Local Lymph Node Assay, LLNA – OECD TG 429), from extreme to non-sensitisers. The adapted method showed full concordance with the Validated Reference Method (VRM) – the KeratinoSens™ published protocol - as follows: Table 1 shows that all 21 chemicals were classified in agreement with the VRM in terms of their skin sensitisation potential. In common with the VRM, the adapted method identified two “false negatives” (phenyl benzoate and eugenol) and one “false positive” (4-methoxyacetophenone) when compared with the LLNA. It is worth noting that the classification of chemicals as “false positive” or “false negative” by *in vitro* tests makes an assumption that the classification generated by the animal *in vivo* test (in this case the LLNA) is predictive of the human response. This is a common practice in the development of *in vitro* methods that is worthy of review.

The Performance Standards stipulate that “The accuracy, sensitivity and specificity of the proposed similar or modified test method should be comparable or better to that of the VRM. The accuracy, sensitivity and specificity obtained with the 20 reference substances should all be equal or higher ( $\geq$ ) than 80.0% (actual for KeratinoSens™ based on the 20 reference substances and using a weighted calculation: 87.0% accuracy, 86.7% sensitivity and 87.5% specificity).”

All acceptance criteria of the Performance Standards were met as follows: Results obtained showed that the accuracy of the adapted method was 85.7%, the sensitivity was 84.6%, and the specificity was 87.5%. The “false positive” and “false negative” results were the same as those obtained with the VRM and no strong or extreme sensitizers were under-predicted. The cell line was the one used in the VRM and the endpoint was the same (luciferase activity measurement).

Our animal-product-free adaptation of the KeratinoSens™ test has proven to be robust and reproducible in our Good Laboratory Practice (GLP) accredited laboratory for over two years, and is performed in full compliance with GLP requirements.

The adapted method provides added value from both scientific and ethical perspectives. From a scientific perspective, the eradication of animal-derived components from *in vitro* test systems, and their replacement with human

equivalents, directly increases the relevance of the test system to humans and would be expected to at least equal, if not enhance, the predictivity of the tests across a wide range of chemical categories. The adaptations also provide cosmetic companies with *in vitro* systems with enhanced ethical value, in response to consumer demands to avoid the use of animals and animal-derived components altogether in cosmetics testing. Human equivalents to animal-derived cell culture components are freely available commercially, and should be sourced from a reputable supplier with full donor consent and safety screening. In our experience, increases in the cost of such products over the animal-derived equivalents are small and do not have an impact on the overall cost of performing the test. We recommend the use of human serum from pooled donors. As with any type of serum, when a new batch is used, an internal validation of the batch including cell morphology, growth rates and  $I_{max} / EC_{1.5}$  values with the positive control and representative reference chemicals should be conducted, with subsequent reservation of successfully performing batches for long term use. Over the past 2.5 years, more than 15 batches of human serum have been used in our laboratory to perform the test, each batch consisting of pooled serum from 60-70 donors. Over this time, the mean induction value of the positive control, cinnamic aldehyde, showed good reproducibility ( $1.84 \pm 0.322$  at a concentration of  $32 \mu\text{M}$ ).

The adaptation to the KeratinoSens™ test method, described here, will contribute to further international harmonisation of hazard and risk assessment. Regulatory acceptance would enable incorporation into the skin sensitisation IATA (Integrated Approach to Testing and Assessment) as a completely animal-product-free *in vitro* test option. We recently obtained clarification from the European Chemicals Agency (ECHA) that data using the adapted method may be used in REACH submissions, provided that the Performance Standards data, demonstrating equivalence with the VRM, is included in the dossier. The path to inclusion in OECD TG442d is already underway. Formal listing in the Test Guideline would mean that the proposed adaptation to the existing method will become globally available for use by test laboratories, for compliance with a range of chemical safety regulations around the world. Subsequent adoption by other laboratories will collectively avoid the unnecessary use of a large volume of animal-derived serum in the build-up to the 2018 REACH deadline and beyond, while maximizing the relevance of the test to human safety.

## References

- DB-ALM (2013). Protocol 155: KeratinoSens™. <https://ecvam-dbalm.jrc.ec.europa.eu/methods-and-protocols/search/155>
- Emter, R., Ellis, G. and Natsch, A. (2010). Performance of a novel keratinocyte-based reporter cell line to screen skin sensitizers *in vitro*. *Toxicol Appl Pharmacol* 245, 281-290. <https://doi.org/10.1016/j.taap.2010.03.009>
- Even, M.S., Sandusky, C.B. and Barnard, N.D. (2006). Serum-free hybridoma culture: ethical, scientific and safety considerations. *Trends Biotechnol* 24, 105-108. <https://doi.org/10.1016/j.tibtech.2006.01.001>
- Jochems, C.E., van der Valk, J.B., Stafleu, F.R. and Baumans, V. (2002). The use of fetal bovine serum: ethical or scientific problem? *Altern Lab Anim* 30, 219-227.
- Lou, H., Du, S., Ji, Q. and Stolz, A. (2006). Induction of AKR1C2 by phase II inducers: identification of a distal consensus antioxidant response element regulated by NRF2. *Mol Pharmacol* 69, 1662-1672. <https://doi.org/10.1124/mol.105.019794>
- Natsch, A., Ryan, C.A., Foertsch, L., et al. (2013). A dataset on 145 chemicals tested in alternative assays for skin sensitization undergoing prevalidation. *J Appl Toxicol* 33, 1337-52. <https://doi.org/10.1002/jat.2868>
- Natsch, A., Bauch, C., Foertsch, L., et al. (2011). The intra- and inter-laboratory reproducibility and predictivity of the KeratinoSens assay to predict skin sensitizers *in vitro*: results of a ring-study in five laboratories. *Toxicol In Vitro* 25, 733-744. <https://doi.org/10.1016/j.tiv.2010.12.014>
- Natsch, A. (2010). The Nrf2-Keap1-ARE toxicity pathway as a cellular sensor for skin sensitizers--functional relevance and a hypothesis on innate reactions to skin sensitizers. *Toxicol Sci* 113, 284-292. <https://doi.org/10.1093/toxsci/kfp228>
- Natsch, A. and Emter, R. (2008). Skin Sensitizers Induce Antioxidant Response Element Dependent Genes: Application to the *In Vitro* Testing of the Sensitization Potential of Chemicals. *Toxicol Sci* 102, 110-119. <http://doi.org/10.1093/toxsci/kfm259>
- OECD (2015). Performance standards for assessment of proposed similar or modified *in vitro* skin sensitization ARE-NRF2 Luciferase test methods. Series on Testing and Assessment No. 213. ENV/JM/MONO(2015)6. <http://bit.ly/2mwKSdA>
- Reisinger K., Hoffmann S., Alépée N., et al. (2015). Systematic evaluation of non-animal test methods for skin sensitisation safety assessment. *Toxicol In Vitro* 29, 259-270. <https://doi.org/10.1016/j.tiv.2014.10.018>

## Acknowledgements

The authors would like to thank Lush Ltd (Poole, UK) for their support which enabled this work.

## Conflict of Interest Statement

None of the authors have potential conflicts of interest.

## Correspondence Address

Dr Carol Treasure  
The Innovation Centre  
Sci-Tech Daresbury  
Keckwick Lane, Daresbury,

Cheshire WA4 4FS, UK  
Phone: +44 (0)1925 607 134.  
Email: [carol.treasure@x-cellr8.com](mailto:carol.treasure@x-cellr8.com)